

Evaluation of a Novel Collagen Dressing Used with a Peel and Place Dressing

Kathleen Derrick M.S., Amy K. McNulty, Ph.D., Kris Kieswetter Ph.D., M.B.A.

Solventum, San Antonio, TX

Introduction

- Negative Pressure Wound Therapy (NPWT) with reticulated open cell foam (ROCF) has advanced wound care by promoting healing.
- Oxidized regenerated cellulose (ORC), a plant-derived material, combined with collagen, forms a dressing that stabilizes growth factors by inhibiting matrix metalloproteinases in the wound environment.

Purpose

- This study assessed the efficacy of a novel collagen/ORC antimicrobial dressing (COA) used with an integrated peel and place NPWT system (PnP) in enhancing wound healing and granulation tissue development in a porcine model.

Methods

- Three domestic swine with full-thickness surgical wounds were treated over seven days with either a standard collagen dressing* or COA.
- Treatments were paired with either traditional NPWT using ROCF dressing[†] with a dressing change on Day 3, or a single-application integrated peel and place dressing[‡].
- The dressing was connected to a negative pressure therapy unit[§] and continuous negative pressure at -125 mmHg was applied to all wounds.
- The study was approved by the Institutional Animal Care and Use Committee (IACUC) and conducted in compliance with applicable regulations.
- At study conclusion, tissue samples were collected for histopathological and morphometric analysis to evaluate granulation tissue quality and thickness.
- Wound volume and area were measured on Days 0, 3, and 7 using 3D photographic imaging.

Results

- No significant differences were found in re-epithelialization percent, bacteria, edema, or serocellular crust thickness between the Collagen/ORC + ROCF with a dressing change at Day 3 and the Collagen/ORC + PnP with a single 7-day dressing placement
- Total wounds treated included: COA + PnP = 6; COA = 4; Collagen/ORC + PnP = 12; Collagen/ORC + ROCF = 6; and PnP = 2.
- Treatments were paired with either traditional NPWT using ROCF[†] with a dressing change on Day 3, or a single-application integrated peel and place dressing[‡].
- 3D photographic imaging for area and volume were not significantly different between all treatment groups (Figures 1-2).




Figure 1. Wound volume (% decrease) from Day 0. Note PnP alone not graphed due to small sample size (n=2).




Figure 2. Wound area (cm²) from Day 0. Note PnP alone not graphed due to small sample size (n=2).

Results (Cont'd)

- Tissues were trimmed per instructions, photographed, processed for paraffin infiltration, embedded, sectioned, mounted on slides, routinely stained with Hematoxylin & Eosin (H&E).
- Granulation tissue maturation was consistent across all treatment groups, with notable collagen deposition for collagen/ORC and COA treated wound sites (Figures 3-5).
- Histopathological findings were similar across treatments, including granulation bed maturity, inflammation, fibrosis, neovascularization, and hemorrhage. Inflammatory markers were not elevated, and no edema or seroma was detected.
- Although major histological differences were not observed, wounds treated with COA + PnP showed a significant difference (p<0.05) in wound area as compared to other treatment groups.
- Wounds treated with PnP had slightly more re-epithelialization than other treatment groups, however, not significant (Figures 6-7).




Figure 3. Peel and Place (PnP). Red line denotes the junction of two combined slides.




Figure 4. L to R: COA, COA + PnP. Red line denotes the junction of two combined slides.




Figure 5. L to R: Collagen/ORC + ROCF; Collagen/ORC + PnP. Red line denotes the junction of two combined slides.

Results (Cont'd)




Figure 6. Wound area (mm) measured by histology (*p<0.05). Note PnP alone not graphed due to small sample size (n=2).




Figure 7. Percent re-epithelialization by histology. Note PnP alone not graphed due to small sample size (n=2).

Conclusions

- No significant differences were found in re-epithelialization percent, bacteria, edema, or serocellular crust thickness between the Collagen/ORC + ROCF with a dressing change at Day 3 and the Collagen/ORC + PnP with a single 7-day dressing placement
- Wounds treated with COA + PnP had significantly smaller histologic area measurement than those in the other treatment groups.
- No adverse effects were observed during the seven-day study period with either collagen dressing (Collagen/ORC and COA) used in conjunction with NPWT at -125 mmHg, whether paired with ROCF and a dressing change or with a single-application PnP dressing.

Presented at the Symposium on Advanced Wound Care Spring/Wound Healing Society, April 8-12, 2026, Charlotte, NC.

NOTE: Specific indications, contraindications, warnings, precautions, and safety information exist for these products and therapies, some of which may be Rx only. Please consult a clinician and product Instructions for Use prior to application.

*Solventum™ Promogran Prisma™ Collagen Matrix with ORC and Silver; †Solventum™ V.A.C.® Therapy with Solventum™ V.A.C.® Granufoam™ Dressing; ‡Solventum™ V.A.C.® Therapy with Solventum™ V.A.C.® Peel and Place Dressing; §Solventum™ ActiV.A.C.™ Therapy Unit (Solventum Corporation, Maplewood, MN); ** inSight®, (eKare Inc., Fairfax, VA).

The authors thank Solventum for poster preparation and production.