



NanoPulse: A Scalable Micromixing Method for Robust and Reproducible RNA-LNP Formulation Across All Scales

Elio Gereige^{a,b}, Claire Counil^a, Jade Giraud^a, Milad Baroud^c, Matthieu Kerhuel^a, Audrey Nsamela^a, Brice Calvignac^b, Chantal Pichon^c, Marie Bonnin^b, Robin Oliveres^a, Thomas Guérinier^a

^a Inside Therapeutics, Bordeaux, France

^b MINT Laboratory – Inserm UMR 1066 and CNRS 6021, University of Angers, Angers, France

^c ART ARNm – Inserm UMS 55 and LI2RSO, University of Orléans, Orléans, France

ABSTRACT

Nanomedicine landscape:

- Growing interest in RNA-LNP therapies, with proven clinical impact.

Translational bottleneck:


- Lack of unified technology to scale LNP production from early-stage R&D to clinical production.
- No single formulation technology ensures identical LNPs across all scales.

NanoPulse technology:

- Only technology that fully covers the whole development pipeline from μL -scale screening to continuous large-scale production.
- First micromixer delivering consistent mixing – thus LNP characteristics, across all scales ($\mu\text{L} \rightarrow \text{L}$) and with different cargos.




TECHNICAL BOTTLENECKS

1. **Small-scale gap:** IJMs and other large-scale mixers struggle to handle discovery volumes ($< 1 \text{ mL}$).
2. **Large-scale gap:** Microfluidics face throughput limits and are prone to aggregation; therefore, are unsuitable for clinical production.
3. **Platform switching & scale-up:** Transitions between methods or platforms alter critical quality attributes (CQAs), requiring re-optimization and re-validation, slowing development.



1. NanoPulse Technology

- **Mixing principle:** Sequential injection of the organic phase (lipids in solvent) and the aqueous phase (with or without cargo, e.g., RNA, DNA) ensures controlled LNP formulation (Figure A). Mixing efficiency validated via semi-Lagrangian CFD simulations (Figure B).
- **Operational control:** Total Flow Rate (TFR), Flow Rate Ratio (FRR), and valve-injection frequency set via the controller (Figure C).
- **Scalability:** Reproducible LNP characteristics from screening to clinical-scale production (goal: 200 μL up to 100+ L).
- **Additional advantages:** Aggregation-free operation and seamless transition from discovery to production without platform switching; reducing regulatory complexity, accelerating process development, minimizing development costs and risks, and ensuring consistent efficacy.



3. Experimental Results

Results in rotating disk format

Efficient formulation (Figure 1):


- mRNA-LNPs visualized at high (167 Hz) and low (12.5 Hz) frequencies using Cryo-TEM.
- Well-formed and intact LNPs with electron-dense internal structure at both frequencies, consistent with successful RNA-LNP assembly.

Continuous operation (Figure 2):

- Stable LNP size & PDI during extended runs, measured using DLS.

Identical RNA-LNP quality at all scales (Figure 3):

- No significant changes in size, PDI, EE%, potency, cell transfection, or protein expression, from 0.25 to 60 mL.



Pulse frequency effect (Figure 4):

- Direct impact on LNP size, aligning with NanoPulse's theoretical model and simulation results.
- Higher injection frequency $\rightarrow$ Smaller LNP size.
- Matching microfluidics (TAMARA) particle sizes by tuning injection frequency.

Adapted to sensitive cargoes (Figure 5):

- Efficient formulation and robust biological performance of saRNA at low and high TFRs.
- Pulse frequency **does not impact** significantly mRNA or saRNA-LNPs potency.
- Protein expression is **comparable to standard microfluidics** (TAMARA) production.




Figure 1

Figure 1

Cryo-TEM of mRNA-LNPs
TFR: 12 mL/min, FRR: 3:1,
[Lipid]: 5 mg/mL, N/P: 6:1
ALC0315:DSPC:Cholesterol:DMG-PEG2000
(50:10:38.5:15)

Figure 2

Size & PDI of empty LNPs
TFR: 36 mL/min, Frequency: 17 Hz,
FRR: 3:1, [Lipid]: 2 mg/mL, V: 4 L
SPC:Cholesterol:DMG-PEG2000
(59:40:1)

Figure 3

siRNA-LNP CQA consistency from 250 μL
to 60 mL
TFR: 20 mL/min, Frequency: 17 Hz,
FRR: 3:1, [Lipid]: 5 mg/mL
ALC0315:DSPC:Cholesterol:DMG-PEG2000 (50:10:38.5:15)

Figure 4

Impact of frequency on siRNA-LNP size
& PDI
TFR: 20 mL/min, FRR: 3:1, [Lipid]: 5 mg/mL
SM102:DSPC:Cholesterol:DMG-PEG200
(50:10:38.5:15)

Figure 5

Effect of RNA modality (saRNA vs
mRNA) and TFR on eGFP protein
expression level in HeLa cells
(Confidential composition)

CONCLUSION

- Controlled LNP formulation via sequential injection of organic and aqueous phases, ensuring efficient mixing, validated by CFD simulations.
- Tunable parameters (TFR, FRR, frequency) for precise LNP size control.
- Seamless scale-up from low-volume research to high-volume continuous production using a single device.
- Reproducible, aggregation-free, volume-independent LNP formulation.

FUTURE DIRECTIONS

- Expand dataset to validate preliminary results.
- Refine the correlation between simulation models and experimental outcomes.
- Characterize CQAs across all scales (morphology, in-vivo performance).
- Combine with downstream processing tools (dilution, TFF, etc.).
- Proceed with a NanoPulse GMP-ready formulation platform.



in collaboration with:



For more info:

contact@insidetx.com



Insidetx.com

