

HIGH-THROUGHPUT ELECTRONEBULIZATION FOR THE SCALABLE PRODUCTION OF CONTROLLED-RELEASE PARTICLES

Ana Salevic¹, Ozan Basar¹, Jose Maria Lagaron^{1,2}

¹Bionanopharma, ²Spanish National Research Council – CSIC, Valencia, Spain
ana.salevic@bionanopharma.com

INTRODUCTION

Polymeric particles:

- Widely used to increase stability and control release of active pharmaceutical ingredients (APIs) in Long-Acting Injectables (LAI) for Subcutaneous (SQ) applications

ELECTRONEBULIZATION¹:

- Proprietary encapsulation technology¹:

Innovative, high-throughput (L/h-scale) nebulization of polymer/API solutions within an electric field at temperatures of 30°C or below

- Main characteristics of the technique:

- controlled particle morphology and architecture, e.g., core-shell
- continuous process
- encapsulated material is collected as a free-flowing powder
- non-thermal process
- ideal for labile APIs, such as biotherapeutics
- scalable to any batch size

METHODS

Particle Production:

- Production of polymeric particles encapsulating various APIs at different throughput rates using the proprietary Capsultek[®] electronebulization technology
- Processing performed at ≤30°C to enable low-temperature particle formation

Particle Formulation:

- Single-polymer and polymer-blend systems
- Various core-shell architectures
- Different processing conditions to tailor particle properties and release profiles

Model Active Pharmaceutical Ingredients (APIs):

- Caffeine
- A commercially available small-molecule drug

Particle Characterization:

- Particle morphology and size distribution by scanning electron microscopy (SEM)

In Vitro Drug Release:


- In vitro* drug-release kinetics evaluated under simulated subcutaneous (SQ) conditions at 37°C
- API release quantified by UV-Visible spectroscopy

References

¹Prieto, C., Evtoski, Z., Pardo-Figueroa, M., Hrakovsky, J., Lagaron, J. M. (2021). Nanostructured valsartan microparticles with enhanced bioavailability produced by high-throughput electrohydrodynamic room-temperature atomization. *Molecular Pharmaceutics*, 18(8), 2947-2958.

RESULTS

- Electronebulization successfully produced submicron- and micron-sized polymeric particles across a broad range of polymer compositions, polymer blends, core-shell configurations, and processing conditions.
- SEM characterization and *In vitro* release studies revealed that both particle morphology and drug-release kinetics were highly dependent on polymer composition, solvent selection, and process parameters.




SQ simulation of caffeine release from the various particle compositions:


Typical particle size distribution by SEM

- PLGA(50:50)/PLGA(50:50): the fastest profile
~70–75% by day 1; plateau at ~86–87% by ~15–20 days
- PLGA(85:15)/PLGA(85:15): more gradual release
~25–30%: day 1, ~75%: ~20 days, ~85%: ~28–30 days
- PLGA(75:25)/PLGA(75:25) and PLLA/PLLA+PEG (5%): intermediate release
increase to ~70–73% by ~28–30 days
- PVA/PLGA(75:25): limited progression after the initial phase
~55–57% over ~5–20 days

Accelerated *in vitro* evaluation of LAI API particles predicting 6-month subcutaneous sustained release:



Burst release suppressed



CONCLUSION

- Electronebulization enabled the high-throughput production of diverse polymeric particles with tunable properties governed by polymer composition and processing conditions.
- By optimizing polymer selection and processing parameters, sustained drug release over extended periods was achieved, demonstrating the potential of this technology as a versatile platform for the development of long-acting injectable (LAI) drug formulations.



BIONANOPHARMA



INSUD PHARMA