

High Throughput Single Particle Characterisation of Lipid Nanoparticles (LNPs) using the Integrated Nanoparticle Size and Payload Analyser (INSPA)




Lan Anh Huynh, Andreas Wallucks, Junjie Bai, and David Juncker
Biomedical Engineering Department,
McGill University, Montréal, Québec, Canada.




Introduction

- Development of **Lipid Nanoparticle (LNP)** formulations for nucleic acid delivery requires laborious optimisation of physical and chemical process parameters to ensure a **high loading efficiency** as well as **high reproducibility of cargo delivery to recipient cells** - an ongoing challenge the LNP field is still actively trying to tackle.¹
- The standard LNP characterisation methods are **DLS** which **lacks single-particle resolution** and is susceptible to bias, and **Cryo-TEM**, which offers single particle resolution **but lacks throughput**.
- Other methods are limited, **i.e. NTA** lacks sensitivity with gradual fall-off for LNPs < 80 nm, while **MRPS** lacks throughput and reliability.
- We introduce the **Integrated Nanoparticle Size and Payload Analyser (INSPA)** that combines **scattering microscopy (iSCAT)** for sizing and **fluorescent staining for RNA loading efficiency quantification**.
- INSPA runs on a conventional fluorescence microscope with a simple beam splitter and uses a microfluidic flow cell to immobilize LNPs on a coverslip for imaging.




1 Current challenge in LNP optimisation



Critical Process Parameters (CPPs): variables during fabrication with significant impact on LNP quality, namely:


- Component molar ratio
- Flowrate Ratio (FRR)
- Total Flowrate (TFR)
- N:P Ratio
- Dialysis/Dilution
- Ionizable lipid
- RNA Concentration
- Aqueous buffer




Large parameter space → optimisation and screening become laborious and can benefit from having a **high throughput characterisation method**.

2 Integrated Nanoparticle Size and Payload Analyzer


An iSCAT-based imaging platform that provides both single particle resolution sizing of LNPs and RNA quantification per LNP.



3 Benchmarking INSPA



4 Introducing a new CQA metric: Per Particle Payload (PPP)



Conclusion & Future work

- INSPA measured both size and loading efficiency of single LNPs down to ~ 30 nm diameter in a few minutes using 15x fold less sample volume than DLS.
- INSPA shortens the design-build-characterize-learn loop, accelerate screening and upon cell transfection help uncover new structure-function relationships thanks to its integrated multiparametric size and payload analysis.
- Future work: Incorporate INSPA into a screening platform for LNP optimisation with a custom fabrication setup, low volume post processing, and low volume cell transfection.

Acknowledgements



References

- McMillan, C. *et al.* Tailoring lipid nanoparticle dimensions through manufacturing processes. *RSC Pharm* **1**, 841-853 (2024). <https://doi.org/10.1039/d4pm00128a>
- Wallucks, A., DeCorwin-Martin, P., Shen, M. L., Ng, A. & Juncker, D. Size photometry and fluorescence imaging of immobilized immersed extracellular vesicles. *Journal of Extracellular Vesicles* **13**, e12512 (2024)
- Ly, H. H., Daniel, S., Soriano, S. K. V., Kis, Z. & Blakney, A. K. Optimization of Lipid Nanoparticles for saRNA Expression and Cellular Activation Using a Design-of-Experiment Approach. *Molecular Pharmaceutics* **19**, 1892-1905 (2022)